



Marina Martínez-Bartolomé, Irene Palacios-Doiztúa, Alberto Garcia-Nieto, Alberto Hernández-Pinto, Rubén Miguez, Antonio Quílez-Álvarez, and Julio Martín\*

Arrays for Cell Nanodevices (A4cell), Madrid, Spain  
\* Presenting author: ✉ [julio.martin@a4cell.com](mailto:julio.martin@a4cell.com)

## SITUATION

- **Live-cell continuous monitoring** of hallmarks for signalling, metabolism, and cell health (e.g. pH, calcium, ROS, oxygen) over long periods of time is priceless.
- Biological samples are **precious and unique**.
- Outputs obtained from the same single cell diminish **variability** and establishes more **reliable correlations and robust statistics**.
- Live-cell imaging in drug discovery and cell biology research is hampered by **drawbacks of traditional fluorescent chemical probes in solution**:
  - cytotoxicity over time,
  - difficulty of multiplexing,
  - metabolization and expelling of molecules from inside the cell,
  - poor cell penetration.

## TARGET


There is a demand for **new technologies advancing high-content live cell imaging assays** from morphological snapshots to physiological feature films.

## PROPOSAL

**Lab-in-a-cell-microchips sensing and filming live cell physiology for drug discovery and cell biology research**




## SPAchip<sup>®</sup> TECHNOLOGY



SPAchip<sup>®</sup>: A silicon device on which multiple highly concentrated fluorescent probes can be covalently printed to provide intracellular readouts over long culture periods


- SPAchip<sup>®</sup> detects pH, Calcium ion and ROS in living single cells over long periods of time, both in **2D and 3D cultures**. *Upcoming O<sub>2</sub> kit.*
- SPAchips are **internalised by cells** after an overnight incubation.
- SPAchip<sup>®</sup> technology is **harmless, non-cytotoxic**.
- **Intracellular and extracellular** on the same sample.
- Amenable to **multiplexing**: 2-sensors-in-one SPAchip.
- Fluorescence intensity is measured from individual SPAchips. **Confocal fluorescence microscopes and HCS/HCA analysers** with at least 20X, yet widefield fluorescence microscopes and flow cytometers as well.

## Cell Imaging Studio





## 3D-Biology: Spheroids

- SPAchip<sup>®</sup> technology enables dynamic, real-time assays monitoring intracellular pH inside 3D cell aggregates over a period of minutes, hours, or even days.
- SPAchips remain in the cytosol reluctant to exocytosis during spheroid handling.



HEK293 spheroids. Timelapse with SpheroCHECK SPAchip<sup>®</sup> pH Single-Detection Kit Green for 24 (A), 28 (B), 48 (C), and 52 hours (D). DRAQ5 in red stained live nuclei.



## Calcium in cell health and receptor signaling



Calcium green SPAchip<sup>®</sup> in HL1 cells responding to the addition of ATP, agonist of purinergic receptors, in the culture media

HL1 cells micrograph at 40X with calcium green detection SPAchip<sup>®</sup> internalized.

## (#3) Precise intra- and extra-cellular pH




**CytoCHECK SPAchip<sup>®</sup> pH RED Single Detection** at different pHs using commercial calibrators. Ratiometric values were obtained by dividing  $\lambda_{emi2}/\lambda_{emi1}$  emission signals in HCS-Operetta with  $\lambda_{exc}=545/15$  nm.


U  
S  
E  
C  
A  
S  
E  
S

## ROS (Reactive Oxygen Species) detection



### OHrad ROS SPAchip<sup>®</sup>



- **KEY FEATURES:**
- **Stability:** active only when processed by intracellular cytosolic esterases
- **Specificity:** hydroxyl radical ( $\cdot\text{OH}$ )
- **Robustness:** pH-insensitive



## Multiplexing pH & Calcium detection



SH-SY5Y cell line (neuroblastoma cells) with nuclei stained in blue and CytoCHECK SPAchip<sup>®</sup> Multi-detection kit in (A) **green channel (Calcium)** and (B) **red channel (pH)**. NB: Both pH and calcium probes are attached to each chip.

## Take Home Messages

- **SPAchip<sup>®</sup> enables real-time monitoring of physiological hallmarks on living single-cells over long periods of time** to gain a deeper understanding of cell biology, cancer and metabolic diseases, and accelerate drug discovery.
- **SPAchip<sup>®</sup> product line comprises pH, Calcium (Single- and Multi-detection kits) and ROS**, and it will be expanded with the **upcoming molecular oxygen kit**.